














VALIDATED DEVICES | COMPATIBILITY














Pac-Dent is dedicated to continuing expansion of the list of devices that are compatible with our Rodin™ 3D Resins. We believe in giving you freedom of choice when it comes to your preferred printing device. See below for our most up-to-date validation and verification status chart.

Validated Status		✓ Verified Compatibility			⌚ Scheduled Validation		● Incompatible	
Compatible Printers		 Rodin® Sculpture 2.0	 Rodin® Titan	 Rodin® Sculpture 1.0	 Rodin® Denture Base 1.0	 Rodin® Denture Base 2.0	 Rodin® Splint 2.0	 Rodin® Clear Bite Splint
	Ackuretta SOL	✓	⌚	✓	✓	✓	✓	✓
	Asiga® MAX (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Asiga® Pico 2 (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Asiga® Pro 4K (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Carbon® 3D (M1, M2, M3)	⌚	⌚	✓	⌚	⌚	⌚	⌚
	DentaFAB SEGA	⌚	✓	✓	✓	⌚	⌚	✓
	Flashforge® Focus 6K	●	✓	✓	✓	✓	✓	✓
	LuxCreo iLux Pro Dental	✓	✓	⌚	⌚	✓	⌚	⌚
	Nexa Xip	⌚	⌚	✓	⌚	✓	✓	⌚
	Phrozen® Sonic Mini 4K	✓	✓	✓	●	✓	✓	●
	Phrozen® Sonic Mini 8K	✓	✓	✓	⌚	✓	✓	✓
	Phrozen® Sonic Mini 8K S	✓	✓	✓	⌚	⌚	⌚	⌚
	Phrozen® Sonic 4K	✓	✓	✓	●	✓	✓	●
	Phrozen® Sonic 4K XL	⌚	⌚	✓	●	✓	✓	●
	RayShape Shape 1+Dental	✓	✓	✓	✓	✓	✓	✓
	RayShape Edge E1	✓	✓	✓	✓	✓	✓	✓
	Shining 3D L4D/L4K	✓	✓	✓	✓	✓	✓	✓
	Shining AccuFab-CEL	⌚	✓	✓	✓	✓	✓	✓
	Whip Mix VeriEKO	⌚	⌚	✓	✓	⌚	⌚	⌚

COMPATIBLE LIGHT CURING DEVICES

Pac-Dent is dedicated to continuing expansion of the list of devices that are compatible with our Rodin™ 3D Resins. We believe in giving you freedom of choice when it comes to your preferred printing device. See below for our most up-to-date validation and verification status chart.

Validated Status		✓ Verified Compatibility			⌚ Scheduled Validation		● Incompatible		
Compatible Light Curing Devices		 Rodin® Sculpture 2.0	 Rodin® Titan	 Rodin® Sculpture 1.0	 Rodin® Denture Base 1.0	 Rodin® Denture Base 2.0	 Rodin® Splint 2.0	 Rodin® Clear Bite Splint	 Rodin® Flexisplint
	Dreve PCU LED N2	✓	✓	✓	✓	✓	✓	⌚	⌚
	Ackuretta Curie Plus	✓	✓	✓	✓	✓	✓	⌚	⌚
	Ackuretta Curie	✓	✓	✓	✓	✓	✓	⌚	⌚
	Formlabs® Form Cure	✓	✓	✓	✓	✓	✓	⌚	⌚
	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓

Compatible Light Curing Devices		 Rodin® Custom Bite Tray	 Rodin® Surgical Guide 2.0	 Rodin® Surgical Guide 1.0	 Rodin® Try-in Denture	 Rodin® Ortho IBT	 Rodin® Model	 Rodin® Castable	 Rodin® Palette
	Dreve PCU LED N2	✓	✓	✓	✓	✓	✓	✓	⌚
	Ackuretta Curie Plus	✓	✓	✓	✓	✓	✓	✓	⌚
	Ackuretta Curie	✓	✓	✓	✓	✓	✓	✓	⌚
	Formlabs® Form Cure	✓	✓	✓	✓	✓	✓	✓	⌚
	Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
FOR BEST AESTHETICS AND MECHANICAL PROPERTIES. MOST RECOMMENDED									