

Abbey Summer Promotions

Take advantage of our summer offers with **discounted prices** and **buy 2 get one free** offers.



Scan to visit our website



Order Online

www.abbeydentalsupplies.com



Order by phone
020 8683 4528



Order on WhatsApp
+44 (0) 7586 208031



Order by Email
sales@abbeydentalsupplies.com

Nu:le Coat Stain & Glaze System

With a rich colour line up and a low viscosity that allows for a thin and uniform abrasion resistant application, the Nu:le Coat Stain and Glaze system allows for the fast, efficient colouring of polymer based crowns.



Nu:le Coat Full Set

Code	Description	Size	Qty	Price
NULEC003	Nu:le Coat Full Set	-	1	£449

Contents Include:

Nu:le Coat Gel	2ml	1
Nu:le Coat Liquid Clear	6ml	1
Nu:le Coat Liquid - White, Grey, Blue, Brown, Yellow, Orange, A Plus, B Plus, C Plus, D Plus	4ml	1
Round Brush	-	1
Flat Brush	-	1
Disposable Plate	-	30pk



Nu:le Coat Gel Single Package

Code	Description	Size	Qty	Price
NULEC017	Nu:le Coat Gel	2ml	1	£58



Nu:le Coat Liquid Single Package

Code	Description	Size	Qty	Price
NULEC018	Nu:le Coat Liquid Clear	6ml	1	£36



Please visit our website or phone our sales team for the full Nu:le Coat colour range.

(All prices exclude VAT and carriage. E & OE)

Tinto Mixing Trays

Tinto mixing trays allow you to store premixed products like stains or paste opaques and keep them protected from dust. Tinto trays space saving unique design provide a small footprint of 118 x 118 x 25 mm. The lid is equipped with a transparent window made out of polycarbonate to allow you to recognize the contents of the tray at a glance.

Two models are available. The first has 9 wells of 18mm diameter and a smaller deeper well for the liquid. The second model has 15 wells of 14mm diameter as well as two small-deep wells for the liquid.

Code	Description	Qty	Price
1300	Tinto, set of 2 trays with 9 and 15 wells	1	£136
1309	Tinto, set of 2 trays with 9 wells	1	£136
1315	Tinto, set of 2 trays with 15 wells	1	£136



(All prices exclude VAT and carriage. E & OE)

RSPCT Brush Kit

RSPCT brushes (100% artificial fibre) have become the next generation of synthetic brushes for thousands of ceramists, who prefer them to natural-hair brushes.

They are now available directly mounted on the new WAVE handles. It's your choice to decide whether you want a set that includes a handle and all 4 head sizes, individual brushes or just handles.

Smile Line RSPCT brushes are used for porcelain build-up, glazing and staining:

- Amazing stability
- Extremely sharp point
- Wonderful resilience
- Perfect water absorption
- Long lasting



#4500-RSPCT-ASST Complete

Code	Description	Price
4500-RSPCT-ASST Complete	<ul style="list-style-type: none">• 1x Complete brush #8Regular,• 1x brush tip of #6Regular, #3, #1	£139

Please visit our website or phone our sales team for the full range of RSPCT replacement tips and handles.

(All prices exclude VAT and carriage. E & OE)

WAVE by Smile Line

WAVE, the new handle from Smile Line with a design that balance's comfort and ergonomics. WAVE handles are compatible with all Smile Line brushes and instruments.

WAVE Handles

Code	Description	Qty	Price
4000-AE-PUR	PUROSANGUE WAVE Handle	1x PUROSANGUE	£23
4000-AE-CAF	CAFÉNOIR WAVE Handle	1x CAFÉNOIR	£23
4000-AE-BON	BONBON ROSE WAVE Handle	1x BONBON ROSE	£23
4000-AE-NAZ	NARAZÉ WAVE Handle	1x NARAZÉ	£23
4000-AE-JAM	JAMES WAVE Handle	1x JAMES	£23

WAVE + RSPCT Kits

Code	Description	Price
4500-RSPCT-ASST-JAM	<ul style="list-style-type: none">1x Complete brush: JAMES Handle + brush tip #8Regular,1x brush tip of #6Regular, #3, #1	£139
4500-RSPCT-ASST-PUR	<ul style="list-style-type: none">1x Complete brush: PUROSANGUE Handle + brush tip #8Regular,1x brush tip of #6Regular, #3, #1	£139
4500-RSPCT-ASST-CAF	<ul style="list-style-type: none">1x Complete brush: CAFÉNOIR Handle + brush tip #8Regular,1x brush tip of #6Regular, #3, #1	£139
4500-RSPCT-ASST-BON	<ul style="list-style-type: none">1x Complete brush: BONBON ROSE Handle + brush tip #8Regular,1x brush tip of #6Regular, #3, #1	£139
4500-RSPCT-ASST-NAZ	<ul style="list-style-type: none">1x Complete brush: NARAZÉ Handle + brush tip #8Regular,1x brush tip of #6Regular, #3, #1	£139

Please visit our website or phone our sales team for the full range of WAVE handles.

(All prices exclude VAT and carriage. E & OE)



#PUROSANGUE



#CAFÉNOIR



#BONBON ROSE



#NAZARÉ



#JAMES



#4500-RSPCT-ASST-PUR

OptiShade Style Italiano

OptiShade is a unique and innovative dental colourimeter. It is specifically designed for measuring the colour of teeth and dental restorations.

OptiShade works exclusively with an iPhone, iPad or iPod (not supplied) via an App available for free on the AppStore.

Set Includes:

- 1 OPTISHADE STYLE ITALIANO unit complete with Capture Guide Cone and Calibration Cap
- 1 lightning cable for connection OPTISHADE STYLE ITALIANO/iPhone-iPad
- 1 USB Charging Cable
- 4 Capture Guide Cones
- 4 Calibration Caps
- 1 Instruction for use booklet and 1 quick start guide OPTISHADE STYLEITALIANO App (free) to be downloaded on Appstore
- 30 days free licence of Matisse software



Code	Description	Qty	Price
70100-OS	OPTISHADE STYLE ITALIANO	1	£1,495 £1,395

(All prices exclude VAT and carriage. E & OE)

Smile Lite MDP2 Master Set

The Smile Lite MDP2 is a plug & play device that is easy to use, doesn't need a special App, calibration process or specific training that will allow you to drastically improve your photographic results with the camera on your smart phone.

- Photos for shade taking
- Every day communication with the lab
- Communication with patients
- Short video clips
- Patients' documentation
- Artistic pictures

It is equipped with three groups of LEDs, each group can be lit up individually. A dimmer function allows you to set four different levels of illumination.



Summer
Special
Price

Free
Blue Tooth
Handle
when ordered
online

Code	Description	Qty	Price
6700-Master	Smile Lite MDP2 Master Set	1	£459

(All prices exclude VAT and carriage. E & OE)



Order Online

www.abbeydentalsupplies.com



Order by phone
020 8683 4528



Order on WhatsApp
+44 (0) 7586 208031



Order by Email
sales@abbeydentalsupplies.com



Order Online

www.abbeydentalsupplies.com



Order by phone
020 8683 4528



Order on WhatsApp
+44 (0) 7586 208031



Order by Email
sales@abbeydentalsupplies.com

Ceramic Fibre Point Polishing Sticks

- Polishing around pit fissures of inlay crowns
- Removal of air bubbles inside the crown or clasp
- Fine adjustment or modification of resin and metal base or attachment
- Shape modification of porcelain



Code	Description	Size (mm)	Qty	Price
YAMAH001	Red Fine	50 x 2.36	1	£15
YAMAH002	Orange Medium	50 x 2.36	1	£15
YAMAH003	Green Coarse	50 x 2.36	1	£15
YAMAH004	Blue Holder for Polishing Sticks	-	1	£5

Dentify Duplicating Gel

SENSE Colourless Duplicating Gel for Gypsum

Agar Agar duplicating gel for the manufacture of plaster models.

- Fast and economical
- Reusable

SENSE Duplicating Gel for Investment and Acrylics

Agar Agar duplicating gel for precise use of investment model and self-cure acrylic.

- High Fitting accuracy
- Reusable



Code	Description	Melting Temp	Working Temp	Colour	Size	Price
640-0306	DENTIFY Sense Duplicating gel for Gypsum	92 - 97°C	48 - 50°C	Colourless	6KG	£80
640-0106	DENTIFY Sense Duplicating gel - for Chrome	90 - 95°C	47 - 48°C	Blue	6KG	£99
640-0206	DENTIFY Sense Duplicating gel - for Chrome	90 - 95°C	47 - 48°C	White	6KG	£77

(All prices exclude VAT and carriage. E & OE)

Shade Peg Original & Shade Peg Release

Shade Peg Original and Release are premium custom peg materials that accurately replicates natural shades of prepared teeth to give a more precise match of your final restoration without the need to remove the crown from the peg and without leaving a hard to remove residue after use.



Code	Description	Colours	Qty	Price
Visit our website or call for codes	Shade Peg Original Syringe	ND1, ND2, ND3, ND4, ND5, ND6, ND7, ND8, ND9, Amalgam Grey, White, Black	10cc	£12
Visit our website or call for codes	Shade Peg Release Syringe	ND1, ND2, ND3, ND4, ND5, ND6, ND7, ND8, ND9, Amalgam Grey, White, Black	10cc	£12
SHADE014	Shade Peg Original Full Kit	Set of 12 tubes	10cc each	£86
SHADE013	Shade Peg Release Full kit	Set of 12 tubes	10cc each	£86

(All prices exclude VAT and carriage. E & OE)

Vanishing Spray

Dental scan spray for extraoral use

SCANTIST 3D – vanishing is a revolutionary scan spray, which was developed for extraoral applications in dental laboratories.

A coating with SCANTIST 3D vanishing spray improves the measurement properties of tooth models, plaster casts, preparations, stumps and impressions.

Product Features

- Self-disappearing, extraoral 3D scan spray for dental applications
- No contamination of object, environment and user
- No cleaning required – cost savings
- Free of titanium dioxide (TiO₂) and other pigments



Code	Description	Size	Qty	Price
250024	Scantist 3D Spray	200ml	1	£52

A·UNO

Direct Composite Resin

Introducing A·UNO, the next stage of direct composite resins. Shade selection is not required as it can match the colour of the tooth by it's unique colour adaptation properties.

A·UNO Single Package

Code	Description	Size	Qty	Price
40800001	A·UNO Universal Basic	4.0g (2mL)	1	£26
40800101	A·UNO Universal St Basic	4.0g (2mL)	1	£26
40800201	A·UNO Low Flow Basic Accessory - Needle Tip: 10 pcs.	2.8g (1.5mL)	1	£26
40800301	A·UNO Low Flow St Basic Accessory - Needle Tip: 10 pcs.	2.8g (1.5mL)	1	£26
40800401	A·UNO Flow Basic Accessory - Needle Tip: 10 pcs.	2.8g (1.5mL)	1	£26
40800501	A·UNO Flow St Basic Accessory - Needle Tip: 10 pcs.	2.8g (1.5mL)	1	£26



#A·UNO Universal Basic



#A·UNO Universal St Basic



#A·UNO Low Flow Basic



#A·UNO Low Flow St Basic



#A·UNO Flow Basic



#A·UNO Flow St Basic

(All prices exclude VAT and carriage. E & OE)

Coming Soon
Pre-order for delivery in July

Special Introductory Price



TMR Z Fill 10.

Direct Composite Resin

The TMR Z Fill 10's unique filler properties allow for both high strength and sustained fluoride release that comes with a comprehensive colour line up.

Coming Soon
Pre-order for delivery in July

TMR Z Fill 10. Single Package

Description	Colours	Size	Qty	Price
TMR-Z Fill 10. Universal	Dentine, A1, A2, A3, A3.5, A4, A5, B1, B2, B3, C2, C3, D2, Opaque A2, Opaque A3, OW, BW, E	3.8g (2mL)	1	£20
TMR-Z Fill 10. Flow	Dentine, A1, A2, A3, A3.5, A4, A5, Opaque A2, Opaque A3, Opaque A3.5, Opaque A4, Opaque A5, OW, BW, E	2.6g (1.5mL)	1	£20
TMR-Z Fill 10. Low Flow	Dentine, A1, A2, A3, A3.5, A4, A5, Opaque A2, Opaque A3, Opaque A3.5, Opaque A4, Opaque A5, OW, BW, E	2.6g (1.5mL)	1	£20
TMR-Z Fill 10. Accessories - Needle tips, TMR Aqua Bond, Applicator brush, disposable plates, Multi primer liquid	-		1	

Please visit our website or phone our sales team for product codes



Special Introductory Price

(All prices exclude VAT and carriage. E & OE)

Making Smiles

Follow Abbey Dental Supplies on Social Media:



Scan to visit our website